

For Research Use Only

# FcZero-rAb® CoraLite® Plus 647 Anti-Human IFN-gamma Rabbit Recombinant Antibody



Catalog Number: CL647-FcA98187

## Basic Information

Catalog Number:

CL647-FcA98187

Size:

100 tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0004

GenBank Accession Number:

BC070256

GeneID (NCBI):

3458

ENSEMBL Gene ID:

ENSG00000111537

UNIPROT ID:

P01579

Full Name:

IFN gamma

Calculated MW:

166 aa, 19 kDa

Purification Method:

Protein A purification

CloneNo.:

241390B10

Recommended Dilutions:

FC (Intra): 5 ul per  $10^6$  cells in a 100  $\mu$ l suspension

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

## Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Positive Controls:

FC (Intra) : PMA and ionomycin treated human PBMCs,

## Background Information

Interferon gamma (IFN  $\gamma$ ), is a type II interferon that provides immunity against bacterial, viral and protozoan infections. In its active form, IFN  $\gamma$  is a glycosylated, non-covalently linked homodimer of 29-32 kDa subunits. It is produced by a number of immune cell types including natural killer cells, natural killer T cells, and effector lymphocyte T cells following antigenic and inflammatory triggers. The IFN  $\gamma$  dimer binds to its cognate receptor which has two subunits: IFN- $\gamma$  R1 which is the ligand-binding chain ( $\alpha$  chain) and IFN- $\gamma$  R2, the signal-transducing chain ( $\beta$  chain). Binding to the receptor activates the JAK/STAT pathway which in turn activates IFN  $\gamma$  responsive genes. While IFN  $\gamma$  can inhibit viral replication, it also works as an immune-modulator and immune-stimulator by increasing surface expression of class I MHC proteins. (PMID: 19268625; 10688427)

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10<sup>6</sup> untreated or PMA and ionomycin treated human PBMCs were intracellularly stained with FcZero-rAb® CoraLite® Plus 488 Anti-Human CD3 (OKT3) and 5 ul CoraLite® Plus 647 Anti-Human IFN-gamma Rabbit RecAb (CL647-FcA98187, Clone: 241390B10). Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).

1x10<sup>6</sup> PMA and ionomycin treated human PBMCs were intracellularly stained with 5 ul CoraLite® Plus 647 Anti-Human IFN-gamma Rabbit RecAb (CL647-FcA98187, Clone: 241390B10) (red) or FcZero-rAb® CoraLite® Plus 647 Rabbit IgG Isotype Control Recombinant Antibody (CL647-FcA98136, Clone: 240953C9) (blue). 1x10<sup>6</sup> untreated human PBMCs were intracellularly stained with 5 ul CoraLite® Plus 647 Anti-Human IFN-gamma Rabbit RecAb (CL647-FcA98187,

