

For Research Use Only

FcZero-rAb® CoraLite® Plus 647 Anti-Human FGFR2/CD332 Rabbit Recombinant Antibody

Catalog Number: CL647-FcA98244



Basic Information

Catalog Number:

CL647-FcA98244

Size:

100 tests, 5 µl/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1265

GenBank Accession Number:

NM_000141.4

GeneID (NCBI):

2263

UNIPROT ID:

P21802

Full Name:

fibroblast growth factor receptor 2

Calculated MW:

92kd

Purification Method:

Protein A purification

CloneNo.:

241383D11

Recommended Dilutions:

FC: 5 µl per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : Transfected HEK-293T cells,

Background Information

FGFR2 (Fibroblast growth factor receptor 2), also known as CD332, is a tyrosine-protein kinase that acts as a cell-surface receptor for fibroblast growth factors and plays an essential role in the regulation of cell proliferation, differentiation, migration, and apoptosis. Ligand binding leads to the activation of several signaling pathways, such as RAS, MAPK1/ERK2, MAPK3/ERK1, and the MAP Kinase signaling pathway, as well as the AKT1 signaling pathway. Mutations in the gene of FGFR2 are associated with Crouzon syndrome, Pfeiffer syndrome, Craniosynostosis, Apert syndrome, Jackson-Weiss syndrome, Beare-Stevenson cutis gyrata syndrome, Saethre-Chotzen syndrome, and syndromic craniosynostosis.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ HEK-293T cells or FGFR2/CD332 transfected HEK-293T cells were surface stained with 5 ul CoraLite® Plus 647 Anti-Human FGFR2/CD332 Rabbit RecAb (CL647-FcA98244, Clone: 241383D11). Cells were not fixed.

1x10⁶ HEK-293T cells or FGFR2/CD332 transfected HEK-293T cells were surface stained with 5 ul CoraLite® Plus 647 Anti-Human FGFR2/CD332 Rabbit RecAb (CL647-FcA98244, Clone: 241383D11) (red) or FcZero-rAb® CoraLite® Plus 647 Rabbit IgG Isotype Control Recombinant Antibody (CL647-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.