

For Research Use Only

CoraLite® Plus 647 Anti-Human CD107a / LAMP1 Rabbit scFv Antibody

Catalog Number: CL647-SF98266-2



Basic Information

Catalog Number: CL647-SF98266-2	GenBank Accession Number: NM_005561.4	Purification Method: Affinity purification
Size: 100ug , Concentration: 500 ug/ml by Nanodrop;	GeneID (NCBI): 3916	CloneNo.: 242116E8
Source: Rabbit	ENSEMBL Gene ID: ENSG00000185896	Recommended Dilutions: IF/ICC: 1:50-1:500 FC (Intra): 0.25 ug per 10 ⁶ cells in a 100 µl suspension
Type: scFv	Full Name: lysosomal-associated membrane protein 1	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Immunogen Catalog Number: EG2611	Calculated MW: 45 kDa	

Applications

Tested Applications: IF/ICC, FC (Intra)	Positive Controls:
Species Specificity: human	IF/ICC : HeLa cells, FC (Intra) : Jurkat cells,

Background Information

LAMP1 (CD107a) is a heavily glycosylated membrane protein enriched in the lysosomal membrane. LAMP1 is extensively glycosylated with asparagine-linked oligosaccharides which protect it from intracellular proteolysis (PMID: 10521503). Although LAMP1 is expressed largely in the endosome-lysosomal membrane of cells, it is also found on the plasma membrane (PMID: 16168398). Elevated LAMP1 expression at the cell surface has also been detected during platelet and granulocytic cell activation, as well as in some tumor cells (PMID: 29085473). LAMP1 functions to provide selectins with carbohydrate ligands. This protein has also been shown to be a marker of degranulation on lymphocytes such as CD8+ and NK cells and may also play a role in tumor cell differentiation and metastasis (PMID: 18835598; 29085473; 9426697).

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
50mM HEPES, 150mM NaCl, 5mM EDTA, 0.1% CHAPS, 10% Sucrose and 0.1% Proclin300, pH7.3

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 Jurkat cells were intracellularly stained with 0.25 μ g CoraLite® Plus 647 Anti-Human CD107a / LAMP1 Rabbit scFv Antibody (CL647-SF98266-2, Clone: 242116E8) (red) or 0.25 μ g CoraLite® Plus 647 Rabbit scFv Isotype Control (CL647-SF98136, Clone: 240953C9) (blue). Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 647 Anti-Human CD107a / LAMP1 Rabbit scFv Antibody (CL647-SF98266-2, Clone: 242116E8) at dilution of 1:200.