

For Research Use Only

CoraLite®647-conjugated CYP1A1 Polyclonal antibody



Catalog Number: CL650-13241

Featured Product

Basic Information

Catalog Number: CL650-13241	GenBank Accession Number: BC023019	Purification Method: Antigen affinity purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 1543	Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Source: Rabbit	Full Name: cytochrome P450, family 1, subfamily A, polypeptide 1	
Isotype: IgG	Calculated MW: 512 aa, 58 kDa	
Immunogen Catalog Number: AG3990	Observed MW: 50-58 kDa	

Applications

Tested Applications:
FC (Intra)

Species Specificity:
human, mouse, pig

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

