

For Research Use Only

CoraLite®647-conjugated p70(S6K) Polyclonal antibody



Catalog Number:CL650-14485

Featured Product

Basic Information

| | | |
|---|---|---|
| Catalog Number: CL650-14485 | GenBank Accession Number: BC053365 | Purification Method: Antigen affinity purification |
| Size: 100ul , Concentration: 1000 µg/ml by Nanodrop; | GeneID (NCBI): 6198 | Excitation/Emission maxima wavelengths: 654 nm / 674 nm |
| Source: Rabbit | Full Name: ribosomal protein S6 kinase, 70kDa, polypeptide 1 | |
| Isotype: IgG | Calculated MW: 59 kDa | |
| Immunogen Catalog Number: AG5883 | Observed MW: 60-75 kDa | |

Applications

Tested Applications:
FC (Intra)

Species Specificity:
human, mouse, rat

Background Information

RPS6KB1(Ribosomal protein S6 kinase beta-1) is also named as STK14A, p70 S6KA and belongs to the S6 kinase subfamily. RPS6KB1 is a major substrate of MTOR and acts as a crucial effector of MTOR signaling pathway. It plays a key role in cell growth and proliferation by regulating INS sensitivity, metabolism, protein synthesis, and cell cycle. RPS6KB1 may play an important role in the progression of HCC and could serve as a potential molecular target for HCC therapy (PMID:22684641). RPS6KB1 is a 70 kDa protein and has 5 isoforms with the calculated molecular mass of 51-59 kDa produced by alternative initiation.

Storage

Storage:
Store at -20°C. Avoid exposure to light.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data