

For Research Use Only

# CoraLite®647 Anti-Human CD11a (TS2/4)



Catalog Number: CL650-65194

## Basic Information

<b>Catalog Number:</b> CL650-65194	<b>GenBank Accession Number:</b> BC008777	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 100tests , 5 µl/test	<b>GeneID (NCBI):</b> 3683	<b>CloneNo.:</b> TS2/4
<b>Source:</b> Mouse	<b>Full Name:</b> integrin, alpha L (antigen CD11A (p180), lymphocyte function-associated antigen 1; alpha polypeptide)	<b>Excitation/Emission maxima wavelengths:</b> 654 nm / 674 nm
<b>Isotype:</b> Mouse IgG1, kappa	<b>Calculated MW:</b> 129 kDa	

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Human

## Background Information

CD11a, also known as integrin αL or LFA-1a, is a 170-180 kDa transmembrane glycoprotein that forms a heterodimer with CD18 (PMID: 7027264; 1672643). CD11a/CD18 (LFA-1) is expressed by all leukocytes and mediates cell adhesion through interactions with its ligands, intercellular adhesion molecule 1 (ICAM-1), ICAM-2, and ICAM-3 (PMID: 7479767). CD11a/CD18 also functions in lymphocyte costimulatory signaling (PMID: 1972160).

## Storage

**Storage:**  
Store at 2-8°C. Avoid exposure to light.

**Storage Buffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data