

For Research Use Only

CoraLite®647 Rat IgG2a Isotype Control (2A3)



Catalog Number:CL650-65209

Basic Information

Catalog Number: CL650-65209	GenBank Accession Number: GeneID (NCBI):	Purification Method: Affinity purification
Size: 100ug, 200 µg/ml	Full Name:	CloneNo.: 2A3
Source: Rat		Excitation/Emission maxima wavelengths: 654 nm / 674 nm
Isotype: IgG2a, κ		

Applications

Tested Applications:
FC
Species Specificity:

Background Information

The 2A3 immunoglobulin is useful as an isotype-matched control. The 2A3 immunoglobulin binds to TNP (trinitrophenol), a hapten not expressed on mouse or human cells. It is used as an isotype control for rat IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:
Store at 2-8°C. Avoid exposure to light.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data