

For Research Use Only

CoraLite®647 Rat IgG2b Isotype Control (LTF-2)



Catalog Number:CL650-65211

Basic Information

Catalog Number:
CL650-65211

Size:
100ug, 200 µg/ml

Source:
Rat

Isotype:
IgG2b

GenBank Accession Number:
GeneID (NCBI):

Full Name:

Purification Method:
Affinity purification

CloneNo.:
LTF-2

Excitation/Emission maxima
wavelengths:
654 nm / 674 nm

Applications

Tested Applications:
FC

Species Specificity:

Background Information

The LTF-2 immunoglobulin is useful as an isotype-matched control. The LTF-2 immunoglobulin has an unknown binding specificity and is used as an isotype control for rat IgG2b antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:
Store at 2-8°C. Avoid exposure to light.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data