

For Research Use Only

CoraLite® Plus 750 Anti-Mouse CD24 (M1/69)



Catalog Number: CL750-65121

Basic Information

Catalog Number: CL750-65121	GenBank Accession Number: BC075622	Purification Method: Affinity purification
Size: 100ug, 0.5 mg/ml	GeneID (NCBI): 12484	CloneNo.: M1/69
Source: Rat	Full Name: CD24a antigen	Excitation/Emission maxima wavelengths: 755 nm / 780 nm
Isotype: IgG2b, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 750 Anti-Mouse CD24 (CL750-65121, Clone:M1/69) (red), or 0.5 ug Isotype Control. Cells were not fixed.