

For Research Use Only

# CoraLite® Plus 750 Mouse IgG2a Isotype Control (C1.18.4)

Catalog Number: CL750-65208



## Basic Information

Catalog Number:

CL750-65208

Size:

100ug, 0.2 mg/mL

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Unknown

CloneNo.:

C1.18.4

Excitation/Emission maxima wavelengths:

755 nm / 780 nm

## Applications

Tested Applications:

FC

Species Specificity:

n/a

## Background Information

The C1.18.4 immunoglobulin is useful as an isotype-matched control. The C1.18.4 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface stained with 0.2 ug CoraLite® Plus 750 Anti-Mouse IgG2a Isotype Control (CL750-65208, Clone:C1.18.4) (red) or CL750 Anti-Human CD3 (OKT3) (green) or unstained (blue). Cells were not fixed. Lymphocytes were gated.