

For Research Use Only

CoraLite® Plus 750 Anti-Mouse CD48 (HM48-1)



Catalog Number: **CL750-65227**

Basic Information

Catalog Number:

CL750-65227

Size:

100ug , 0.5 mg/ml

Source:

Armenian Hamster

Isotype:

IgG

GenBank Accession Number:

BC060977

GeneID (NCBI):

12506

UNIPROT ID:

P18181

Full Name:

CD48 antigen

Purification Method:

Affinity purification

CloneNo.:

HM48-1

Excitation/Emission maxima
wavelengths:

755 nm / 780 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 750 Anti-Mouse CD48 (CL750-65227, Clone: HM48-1) (red), or 0.5 ug Control Antibody. Cells were not fixed.