

For Research Use Only

CoraLite® Plus 750 Anti-Rat CD43 (W3/13)



Catalog Number: **CL750-65237**

Basic Information

Catalog Number: CL750-65237	GenBank Accession Number: GeneID (NCBI): 24796	Purification Method: Affinity purification
Size: 100ug, 0.5 mg/ml	Full Name: sialophorin	CloneNo.: W3/13
Source: Mouse	Calculated MW: 38 kDa	Excitation/Emission maxima wavelengths: 755 nm / 780 nm
Isotype: IgG1, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Rat

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ rat splenocytes cells were surface stained with 0.5 ug CoraLite® Plus 750 Anti-Rat CD43 (CL750-65237, Clone:W3/13), or 0.5 ug CoraLite® Plus 750 Mouse IgG1 Isotype Control (MOPC-21) (CL750-65124, Clone: MOPC-21). Cells were not fixed.