

For Research Use Only

CoraLite® Plus 750 Mouse IgM Isotype Control (MOPC-104E)



Catalog Number: CL750-65243

Basic Information

Catalog Number: CL750-65243	GenBank Accession Number: GeneID (NCBI):	Purification Method: N/A
Size: 100ug, 200 µg/ml	Full Name:	CloneNo.: MOPC-104E
Source: Mouse		Excitation/Emission maxima wavelengths: 755 nm / 780 nm
Isotype: IgM, lambda		

Applications

Tested Applications:
FC
Species Specificity:

Background Information

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human whole blood cells were surface stained with 0.2 ug CoraLite® Plus 750 Anti-Mouse IgM Isotype Control (red), 5 ul CoraLite 750 Plus Anti-Human STRO-1 (green), or unstained (blue). Cells were not fixed. CD45dim, CD34+ lymphocytes were gated.