

For Research Use Only

# CoraLite® Plus 750 Mouse IgG2a Isotype Control (MOPC-173)



Catalog Number: CL750-65244

## Basic Information

Catalog Number:  
CL750-65244

Size:  
100ug, 0.2 mg/ml

Source:  
Mouse

Isotype:  
IgG2a, kappa

GenBank Accession Number:  
GeneID (NCBI):

Full Name:

Purification Method:  
Affinity purification

CloneNo.:  
MOPC-173

Excitation/Emission maxima  
wavelengths:  
755 nm / 780 nm

## Applications

Tested Applications:  
FC

Species Specificity:

## Background Information

### Storage

Storage:  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface stained with 0.2 ug CoraLite® Plus 750 Anti-Mouse IgG2a Isotype Control (CL750-65244, Clone:MOPC-173) (red), Anti-Human CD3 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.