

For Research Use Only

CoraLite® Plus 750 Anti-Mouse GL7 (GL-7)



Catalog Number: **CL750-65261**

Basic Information

Catalog Number:
CL750-65261

Size:
100ug, 0.5 mg/ml

Source:
Rat

Isotype:
IgM, kappa

GenBank Accession Number:

GeneID (NCBI):
104230

Full Name:
lymphocyte antigen 77

Purification Method:

Affinity purification

CloneNo.:

GL-7

Excitation/Emission maxima
wavelengths:

755 nm / 780 nm

Applications

Tested Applications:
FC

Species Specificity:
Human, Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.


Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ unstimulated or anti-CD3/CD28 stimulated mouse splenocytes were surface co-stained with CoraLite® Plus 488 Anti-Mouse CD3 (17A2) and 0.5 ug CoraLite® Plus 750 Anti-Mouse GL7 (CL750-65261, Clone:GL-7). Cells were not fixed.