

For Research Use Only

CoraLite® Plus 750 Mouse IgG2a Isotype Control (MOPC-173) Recombinant Antibody

Catalog Number: CL750-65571



Basic Information

Catalog Number:

CL750-65571

Size:

100ug , 500 ug/ml

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

MOPC-173

Excitation/Emission maxima
wavelengths:

755 nm / 780 nm

Applications

Tested Applications:

FC

Species Specificity:

n/a

Positive Controls:

FC : human PBMCs,

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



1×10^6 human PBMCs were surface stained with 0.2 μ g CoraLite® Plus 750 Mouse IgG2a Isotype Control (MOPC-173) RecAb (CL750-65571, Clone: MOPC-173) (red), or CoraLite® Plus 750 Anti-Human CD3 (OKT3) (CL750-65133, Clone: OKT3) (green), or unstained (blue). Cells were not fixed.