

For Research Use Only

CoraLite® Plus 750-conjugated MAPKAPK2 Recombinant antibody

Catalog Number: CL750-82929-2



Basic Information

Catalog Number:

CL750-82929-2

Size:

100ul , Concentration: 1000 ug/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34196

GenBank Accession Number:

BC036060

GeneID (NCBI):

9261

UNIPROT ID:

P49137

Full Name:

mitogen-activated protein kinase-
activated protein kinase 2

Calculated MW:

400 aa, 46 kDa

Observed MW:

42-46 kDa

Purification Method:

Protein A purification

CloneNo.:

230197D6

Recommended Dilutions:

WB 1:500-1:1000

Excitation/Emission maxima
wavelengths:

755 nm / 780 nm

Applications

Tested Applications:

WB

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, HeLa cells

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with CL750-82929-2 (MAPKAPK2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.