

For Research Use Only

# FcZero-rAb® CoraLite® Plus 750 Anti-Mouse CD90.2 (30-H12) Rabbit Recombinant Antibody

Catalog Number: CL750-FcA65656



## Basic Information

**Catalog Number:**  
CL750-FcA65656

**Size:**  
100ug, 500 ug/ml

**Source:**  
Rabbit

**Isotype:**  
IgG

**GenBank Accession Number:**  
BC054436

**GeneID (NCBI):**  
21838

**UNIPROT ID:**  
P01831

**Full Name:**  
thymus cell antigen 1, theta

**Purification Method:**  
Protein A purification

**CloneNo.:**  
30-H12

**Recommended Dilutions:**  
FC: 0.25 ug per 10<sup>6</sup> cells in a 100 µl suspension

**Excitation/Emission maxima wavelengths:**  
755 nm / 780 nm

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
mouse

**Positive Controls:**

FC : mouse thymocytes,

## Background Information

CD90 (Thy-1) is a 25 kDa, GPI-linked membrane glycoprotein that belongs to immunoglobulin superfamily (PMID: 6177036; 6153212). Originally described as a brain thymus cross-reactive antigen, it is found in large quantities on mouse and rat thymocytes and central nervous system cells (PMID: 83175). CD90 has been postulated to be involved in cellular recognition, adherence, and T cell activation (PMID: 7683034). Mouse CD90 (Thy-1) is defined as a differentiation alloantigen, represented as two serologically distinguishable allelic forms, designated CD90.1 (Thy-1.1) and CD90.2 (Thy-1.2) (PMID: 83175). CD90.2 is expressed by thymocytes, T cells, neurons and hematopoietic stem cells in most strains of mice (C3H, BALB/c, C57BL/6, DBA, and others), whereas CD90.1 is expressed in some mouse strains like AKR (PMID: 5919593).

## Storage

**Storage:**

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1x10<sup>6</sup> mouse thymocytes were surface stained with 0.25 ug CoraLite® Plus 750 Anti-Mouse CD90.2 (30-H12) Rabbit RecAb (CL750-FcA65656, Clone: 30-H12) or 0.25 ug Rabbit IgG Isotype Control RecAb (CL750-FcA98136, Clone: 240953C9). Cells were not fixed.



1x10<sup>6</sup> mouse thymocytes were surface stained with 0.25 ug CoraLite® Plus 750 Anti-Mouse CD90.2 (30-H12) Rabbit RecAb (CL750-FcA65656, Clone: 30-H12) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (CL750-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.