

For Research Use Only.
Not For Use In Diagnostics.

Streptavidin-PerCP Conjugate

Catalog Number: CP-PF00030



Information

Catalog Number:

CP-PF00030

Host:

Applications:

FC, IF

Reactivity:

Physical State:

Liquid

Conjugation:

PerCP

Recommended Dilutions

1:200-1:1000 for flow cytometry and immunofluorescence

Storage

Storage:

2-8°C for one year.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.4.

Aliquoting is unnecessary for -20°C storage

Purity


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Flow cytometric analysis of human PBMCs using Biotin anti-human CD8 and PerCP Streptavidin. Human PBMCs were stained with Biotin CD8 Recombinant Antibody (Cat.No. Biotin-65559, Clone: UCHT4), washed with PBS, and stained with PerCP Streptavidin (Cat.No. CP-PF00030) and analyzed by flow cytometry.