

For Research Use Only.
Not For Use In Diagnostics.

PerCP-Cyanine5.5-conjugated Streptavidin Rabbit RecAb

Catalog Number: CPY5-PF00030



www.ptglab.com

Information

Catalog Number:

CPY5-PF00030

Host:

Applications:

FC, IF

Reactivity:

Physical State:

Liquid

Conjugation:

PerCP-Cyanine5.5

Recommended Dilutions

1:200-1:1000 for flow cytometry and immunofluorescence

Storage

Storage:

2-8°C for one year.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.4.

Aliquoting is unnecessary

Purity


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB
(1-888-478-4522)
(toll free in USA),
or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with Biotin Anti-Human CD8 (UCHT4) Mouse IgG2a Recombinant Antibody (Biotin-65559, Clone: UCHT4), washed with PBS, and stained with PerCP-Cyanine5.5 Streptavidin Recombinant Antibody (CPY5-PF00030) and analyzed by flow cytometry. Cells were not fixed. Lymphocytes were gated.