

For Research Use Only

Recombinant Human IGF1 protein (hFc Tag)



Catalog Number: Eg0206

Basic Information

Species:
Human

Purity:
>90 %, SDS-PAGE

Tag:
hFc Tag

EC50:
20-80 ng/mL

Technical Specifications

Purity:
>90 %, SDS-PAGE

Endotoxin Level:
<0.1 EU/μg protein, LAL method

Source:
HEK293-derived Human IGF1 protein Gly49-Ala118 (Accession# P05019-1) with a human IgG1 Fc tag at the N-terminus.

GeneID:
3479

Accession:
P05019-1

Predicted Molecular Mass:
33.7 kDa

SDS-PAGE:
34-36 kDa, reducing (R) conditions

Formulation:
Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Immobilized Human IGF1 R (GST tag) at 1 μg/mL (100 μL/well) can bind Human IGF1 (hFc tag) with a linear range of 20-80 ng/mL.

Storage and Shipping

Storage:
It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:
The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

Insulin-like-growth factor 1 (IGF1), a 70 amino-acid peptide hormone is the principal mediator of biochemical effects of growth hormone (GH). IGF1 is an important growth factor in the regulation of cell proliferation and differentiation. IGF1 is largely synthesized in the liver (75%) and, to a lesser extent, in peripheral tissues. IGF1 have been shown to play an essential role in preventing the formation of fatty liver. IGF1 is a potent mitogen and is inhibited by IGF-binding protein-3 (IGFBP3). High serum IGF1 and low IGFBP3 are associated with increased risk of several carcinomas.

References

1. Adamo ML. et al. (1993) Adv Exp Med Biol. 343:1-11.
2. Adachi Y. et al. (2019) J Gastroenterol Hepatol. 2019 Jun 3.
3. Anisimov VN. et al. (2019) Crit Rev Oncol Hematol. 87(3):201-23.
4. Sonntag WE. et al. (2012) J Gerontol A Biol Sci Med Sci. 67(6):587-98.

Synonyms


IGF1, IGF-1, IGF I, IGFI, IGF-I

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Human IGF 1 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) conditions and stained using Coomassie blue.



Immobilized Human IGF1 R (GST tag) at 1 µg/mL (100 µL/well) can bind Human IGF1 (hFc tag) with a linear range of 20-80 ng/mL