

For Research Use Only


Recombinant Mouse LIF protein (His Tag)



Catalog Number: Eg0732

Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: His Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/µg protein, LAL method Source: HEK293-derived Mouse LIF protein Ser24-Phe203 (Accession# P09056) with a His tag at the C-terminus. GeneID: 16878 Accession: P09056 Predicted Molecular Mass: 20.9 kDa SDS-PAGE: 30-45 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Leukemia inhibitory factor (LIF) is a multi-functional cytokine of the interleukin-6 (IL-6) superfamily. It utilises a receptor that consists of the LIF receptor B and gp130 and this receptor complex is also used by ciliary neurotrophic growth factor (CNTF), oncostatin M, cardiotrophin1 (CT1) and cardiotrophin-like cytokine (CLC). Functionally, LIF has been implicated in a many physiological processes including the induction of hematopoietic differentiation in normal and myeloid leukemia cells, induction of neuronal cell differentiation, regulator of mesenchymal to epithelial conversion during kidney development, and may also have a role in immune tolerance at the maternal-fetal interface.		
References	1. D J Hilton.et al. (1991) J Cell Biochem.46(1):21-6. 2. Nicos A.et al. (2005) Cytokine & Growth Factor Reviews.26(5):533-544. 3. C J Auernhammer.et al. (2000) Endocr Rev.21(3):313-45.		
Synonyms	Lif, leukemia inhibitory factor		
For technical support and original validation data for this product please contact			
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)		E: proteintech@ptglab.com W: ptglab.com	
		This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.	

Selected Validation Data



Purity of Recombinant Mouse LIF was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) conditions and stained using Coomassie blue.