

For Research Use Only

Recombinant Human BCMA/TNFRSF17 protein (hFc Tag)



Catalog Number: Eg0812

Basic Information

Species:
Human

Purity:
>95 %, SDS-PAGE

Tag:
hFc Tag

EC50:
2-8 ng/mL

Technical Specifications

Purity:
>95 %, SDS-PAGE

Endotoxin Level:
<0.1 EU/μg protein, LAL method

Source:
HEK293-derived Human BCMA protein Met1-Ala54 (Accession# Q02223) with a human IgG1 Fc tag at the C-terminus.

GeneID:
608

Accession:
Q02223

Predicted Molecular Mass:
32.4 kDa

SDS-PAGE:
35-40 kDa, reducing (R) conditions

Formulation:
Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Immobilized Human BAFF (GST tag) at 1 μg/mL (100 μL/well) can bind Human BCMA (hFc tag) with a linear range of 2-8 ng/mL.

Storage and Shipping

Storage:
It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:
The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

BCMA (B cell maturation antigen), also known as TNFRSF17, is 20.2-kDa type III transmembrane glycoprotein and is a member of the TNF-receptor superfamily. This receptor is preferentially expressed in mature B lymphocytes and plasma cells, which may be important for B cell development and autoimmune response. BCMA has two agonist ligands: a proliferation-inducing ligand (APRIL) and BAFF. When BCMA binds to APRIL, it transmits signals of cell survival and proliferation; when BCMA binds to BAFF, it mediates the activation of NF-κappaB and MAPK8/JNK. It has been found that the overexpression and activation of BCMA are associated with multiple myeloma (MM) in preclinical models and humans, supporting its potential utility as a therapeutic target for MM.

References

1. Bo Yu et al (2020). J Hematol Oncol. Sep 17;13(1):125.
2. C Madry et al (1998). Int Immunol. Nov;10(11):1693-702.
3. Yu-Tzu Tai et al (2016). Blood. Jun 23;127(25):3225-36.
4. Fabrice Jardin (2022). Biomedicines. Sep 1;10(9):2153.
5. Nina Shah et al (2020). Leukemia. Apr;34(4):985-1005.

Synonyms


BCMA/TNFRSF17, TNFRSF17, B-cell maturation protein, BCM, BCMA

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Human BCMA was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



Immobilized Human BAFF (GST tag) at 1 µg/mL (100 µL/well) can bind Human BCMA (hFc tag) with a linear range of 2-8 ng/mL.