

For Research Use Only

Recombinant Mouse CD200 protein (hFc Tag)(HPLC verified)



Catalog Number: Eg0919

Basic Information

Species:
Mouse

Purity:
>95 %, SDS-PAGE
SEC-HPLC

Tag:
hFc Tag
EC50:
8-32 ng/mL

Technical Specifications

Purity:
>95 %, SDS-PAGE
>90 %, SEC-HPLC

Endotoxin Level:
<0.1 EU/μg protein, LAL method

Source:
HEK293-derived Mouse CD200 protein Gln31-Gly232 (Accession# O54901) with a human IgG1 Fc tag at the C-terminus.

GeneID:
17470

Accession:
O54901

Predicted Molecular Mass:
48.5 kDa

SDS-PAGE:
60-68 kDa, reducing (R) conditions

Formulation:
Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Immobilized Mouse CD200R1 (His tag) at 1 μg/mL (100 μL/well) can bind Mouse CD200 (hFc tag) with a linear range of 8-32 ng/mL.

Storage and Shipping

Storage:
It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:
The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

CD200, also known as OX2, is a type I transmembrane glycoprotein that belongs to the immunoglobulin superfamily. It contains two extracellular immunoglobulin domains, a transmembrane and a cytoplasmic domain. CD200 is expressed on a variety of cell types, including thymocytes, B lymphocytes, a subset of T lymphocytes, follicular dendritic cells, neurons, and endothelial cells. CD200 binds to its receptor, CD200R, which is primarily expressed by myeloid and T cell lineage. The CD200-CD200R interaction plays a role in immunosuppression and suppression of anti-tumor immune responses.

References

- 1.G W McCaughan, et al. (1987) Immunogenetics. 25(5):329-35.
- 2.G J Wright, et al. (2000) Immunity. 13(2):233-42.
- 3.Reginald Gorczynski, et al. (2004) J Immunol. 172(12):7744-9.
- 4.T P Rygiel, et al. (2012) Oncogene. 31(24):2979-88.

Synonyms


CD200 antigen, Mox2, OX2

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Purity of Recombinant Mouse CD200 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



Immobilized Mouse CD200R1 (His tag) at 1 µg/mL (100 µL/well) can bind Mouse CD200 (hFc tag) with a linear range of 8-32 ng/mL.



The purity of Mouse CD200 was greater than 90% as determined by SEC-HPLC.