

For Research Use Only

Recombinant Rabbit CD4 protein (His Tag)



Catalog Number: Eg0927

Basic Information	Species: Rabbit	Purity: >90 %, SDS-PAGE	Tag: His Tag
-------------------	--------------------	----------------------------	-----------------

Technical Specifications

Purity:
>90 %, SDS-PAGE

Endotoxin Level:
<0.1 EU/μg protein, LAL method

Source:
HEK293-derived Rabbit CD4 protein Lys26-Pro396 (Accession# P46630) with a His tag at the C-terminus.

GeneID:
100009152

Accession:
P46630

Predicted Molecular Mass:
42.1 kDa

SDS-PAGE:
39-45 kDa, reducing (R) conditions

Formulation:
Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:
It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:
The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

CD4 is a transmembrane glycoprotein of the immunoglobulin super-family (IgSF), containing four extracellular Ig-like domains, a hydrophobic transmembrane region and a 40-residue cytoplasmic tail. CD4 is expressed T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells. CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction and it is involved in T cell differentiation, development and activation by interacting with antigen-presenting cells (APCs). The cytoplasmic tail of CD4 associates with the Lck kinase, which in turn activates the signaling components of the T cell receptor (TCR)-CD3 complexes. CD4 also serves as a receptor for human immunodeficiency virus (HIV) by binding to the viral envelope glycoprotein gp120.

References

1. Tseng WY, et al. (2013). PLoS One. 8(5):e63963.
2. Bowers K, et al. (1997). Int J Biochem Cell Biol. 29(6):871-875.
3. Graziani-Bowering G, et al. (2002). Exp Cell Res. 279(1):141-152.
4. Moore JP, et al. (1997). Curr Opin Immunol. 9: 551-562.

Synonyms


CD4, T-cell surface antigen T4/Leu-3, T-cell surface glycoprotein CD4

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Rabbit CD4 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.