

For Research Use Only

Recombinant Mouse CD83 protein (rFc Tag)



Catalog Number: Eg1387

Basic Information

Species:
Mouse

Purity:
>90 %, SDS-PAGE

Tag:
rFc Tag

Technical Specifications

Purity:

>90 %, SDS-PAGE

Endotoxin Level:

<0.1 EU/μg protein, LAL method

Source:

HEK293-derived Mouse CD83 protein Met22-Arg133 (Accession# O88324) with a rabbit IgG Fc tag at the C-terminus.

GeneID:

12522

Accession:

O88324

Predicted Molecular Mass:

38.6 kDa

SDS-PAGE:

42-55 kDa, reducing (R) conditions

Formulation:

Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

CD83 is a member of the immunoglobulin (Ig) superfamily, consisting of an extracellular V-type Ig-like domain, a transmembrane domain, and a cytoplasmic tail. CD83 is one of the most prominent markers for fully matured dendritic cells (DCs). It is also found on the surface of other activated hematopoietic cells, including lymphocytes, monocytes, macrophages, and neutrophils. CD83 regulates the maturation, activation, and homeostasis of numerous immune cells. CD83 can also be expressed as a soluble form. Soluble CD83 can bind to DCs and inhibit their maturation.

References

1. Prectel AT. et al. (2007). Arch Dermatol Res. 299(2):59-69.
2. Li Z. et al. (2019). Front Immunol. 10:1312.
3. Grosche L. et al. (2020). Front Immunol. 11:721.
4. Lechmann M. et al. (2002). Int Arch Allergy Immunol. 129(2):113-8.

Synonyms


CD83 antigen, mCD83

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Mouse CD83 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.