

For Research Use Only


Recombinant Mouse CD98/SLC3A2 protein (rFc Tag)(HPLC verified)

Catalog Number: Eg1547




Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE >90 %, SEC-HPLC	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE >90 %, SEC-HPLC Endotoxin Level: <0.1 EU/µg protein, LAL method Source: HEK293-derived Mouse CD98 protein Ala100-Ala526 (Accession# P10852-1) with a rabbit IgG Fc tag at the N-terminus. GeneID: 17254 Accession: P10852-1 Predicted Molecular Mass: 74.8 kDa SDS-PAGE: 75-100 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	CD98 (SLC3A2) is a solute carrier family protein that is a component of the amino acid transporter and pairs with another solute carrier family protein, such as L-type amino acid transporter 1 (LAT1/SLC7A5) or cystine/glutamate transporter (xCT/SLC7A11). CD98 is known to play a role in delivering amino acids into and out of cells as an amino acid transporter. CD98 is highly upregulated in various cancer cells and promotes cancer cell survival, migration, invasion, angiogenesis and drug resistance.		
References	1. Lee Y. et al. (2019). Nat Struct Mol Biol. 26:510–517. 2. Xia P. et al. (2023). Front Oncol. 13:1251100. 3. Lim K. et al. (2024). Cancer Immunol Immunother. 73(11):231.		
Synonyms	4F2 cell-surface antigen heavy chain, 4F2hc, Amino acid transporter heavy chain SLC3A2, CD98, Mdu1		
For technical support and original validation data for this product please contact			
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)		E: proteintech@ptglab.com W: ptglab.com	
This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.			

Selected Validation Data



Purity of Recombinant Mouse CD98 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Mouse CD98/SLC3A2 was greater than 90% as determined by SEC-HPLC.