

For Research Use Only

Recombinant Mouse BAFF protein (rFc Tag)



Catalog Number: Eg1555

Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	<p>Purity: >90 %, SDS-PAGE</p> <p>Endotoxin Level: <1.0 EU/μg protein, LAL method</p> <p>Source: HEK293-derived Mouse BAFF protein Ala127-Leu309 (Accession# Q9WU72) with a rabbit IgG Fc tag at the N-terminus.</p> <p>GeneID: 24099</p> <p>Accession: Q9WU72</p> <p>Predicted Molecular Mass: 47.8 kDa</p> <p>SDS-PAGE:</p> <p>Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.</p>		
Biological Activity	Not tested		
Storage and Shipping	<p>Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.</p> <ul style="list-style-type: none">• Until expiry date, -20°C to -80°C as lyophilized proteins.• 3 months, -20°C to -80°C under sterile conditions after reconstitution. <p>Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.</p>		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background			
References			
Synonyms	B-cell-activating factor, CD257, Tnfsf13b, Tumor necrosis factor ligand superfamily member 13B, Tumor necrosis factor ligand superfamily member 13b, membrane form		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

