

For Research Use Only

Recombinant Human CD72 protein (rFc Tag)



Catalog Number: Eg1755


Basic Information	Species: Human	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Human CD72 protein Arg117-Asp359 (Accession# P21854) with a rabbit IgG Fc tag at the C-terminus. GeneID: 971 Accession: P21854 Predicted Molecular Mass: 54.4 kDa SDS-PAGE: 55-65 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	CD72 is a type II transmembrane glycoprotein of 39-45 kDa and is expressed as a disulphide-linked homodimer. CD72 is a pan B-cell marker covering the full range of differentiation from the very early stages of B cells to mature surface Ig⁺ B cells. It is expressed in all stages of B cell development except plasma cells. CD72 is also expressed by dendritic cells, tissue macrophages in the red pulp of the spleen and von Kupffer cells in the liver. CD72 is the ligand for CD5. CD72 is a negative regulator of BCR signaling. Interaction of CD100 with CD72 strictly tunes the strength of BCR signals and is essential for maintaining immunological homeostasis as well as generating a proper immune response.		
References	1. H Van de Velde, et al. (1991) Nature. 351(6328):662-5. 2. R Schwarting, et al. (1992) Am J Hematol. 41(3):151-8. 3. G Massard, et al. (1996) Ann Thorac Surg. 61(1):252-8. 4. T Adachi, et al. (1998) J Immunol. 160(10):4662-5. 5. Atsushi Kumanogoh, et al. (2005) Int Immunol. 17(10):1277-82.		
Synonyms	B-cell differentiation antigen CD72, CD72 molecule, CD72b, Lyb 2, LYB2		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Human CD72 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.