

For Research Use Only

Recombinant Human Beta-2-Microglobulin protein (rFc Tag)



Catalog Number: Eg2283


Basic Information	Species: Human	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Human Beta-2-Microglobulin protein Ile21-Met119 (Accession# P61769) with a rabbit IgG Fc tag at the C-terminus. GeneID: 567 Accession: P61769 Predicted Molecular Mass: 38.0 kDa SDS-PAGE: 37-45 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Beta-2-microglobulin (B2M) is a 12-kDa protein comprised of 99 amino acids. B2M is a component of MHC class I molecules, which are present on the surface of nearly all nucleated cells. It can be found in body fluids under physiologic conditions as a result of shedding from cell surfaces or intracellular release. B2M has various biological functions, including antigen presentation. Serum B2M is regarded as a marker of disease severity in renal injury, infections, amyloidosis, and aging-related diseases. Investigations reveal that increased synthesis and release of B2M are present in several malignant diseases including multiple myeloma, lymphoma, and solid tumors.		
References	1. D Güssow, et al. (1987) J Immunol. 139(9):3132-8. 2. Jin Xie, et al. (2003) Blood. 101(10):4005-12. 3. Takeo Nomura, et al. (2014) Anticancer Agents Med Chem. 14(3):343-52.. 4. Hanbing Wang, et al. (2021) Cancer Lett. 517:96-104.		
Synonyms	B2M, beta 2 microglobulin, beta 2-Microglobulin, Beta-2-microglobulin form pl 5.3, CDABP0092		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Human Beta-2-Microglobulin was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.