

For Research Use Only

Recombinant Mouse CD302 protein (rFc Tag) (HPLC verified)



Catalog Number: Eg2660

Basic Information

Species:
Mouse

Purity:
>90 %, SDS-PAGE
>90 %, SEC-HPLC

Tag:
rFc Tag

Technical Specifications

Purity:
>90 %, SDS-PAGE
>90 %, SEC-HPLC

Endotoxin Level:
<0.1 EU/μg protein, LAL method

Source:
HEK293-derived Human CD302 protein Asp21-His156 (Accession# Q9DCG2-2) with a rabbit IgG Fc tag at the C-terminus.

GeneID:
66205

Accession:
Q9DCG2-2

Predicted Molecular Mass:
41.6 kDa

SDS-PAGE:
40-47 kDa, reducing (R) conditions

Formulation:
Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

CD302 is a C-type lectin receptor involved in cell adhesion and migration, as well as endocytosis and phagocytosis. In mice and humans, CD302 is primarily expressed in the liver, lungs, lymph nodes (LNs), spleen, and bone marrow. In CD302 gene knockout (CD302KO) mice, the migration frequency and number of myeloid dendritic cells (mDC) in CD302KO lymph nodes were reduced compared to the wild type, confirming the functional role of CD302 in mDC migration.

References

- 1.Lo TH, et al. (2016) J Immunol.197(3):885-98.
- 2.Kato M, et al. (2007) J Immunol. 179(9):6052-63.

Synonyms


CD302 antigen, Clec13a, C-type lectin domain family 13 member A, Type I transmembrane C-type lectin receptor DCL-1

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Mouse CD302 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Mouse CD302 was greater than 90% as determined by SEC-HPLC.