

For Research Use Only

Recombinant Mouse CD37 protein (rFc Tag)



Catalog Number: Eg2729


Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Mouse CD37 protein Arg112-Asn241 (Accession# Q61470) with a rabbit IgG Fc tag at the C-terminus. GeneID: 12493 Accession: Q61470 Predicted Molecular Mass: 41.2 kDa SDS-PAGE: 46-60 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	CD37 (TSPAN26) is a cell surface glycoprotein that is a member of the transmembrane tetrasuperfamily of proteins (TM4SF), which are predominantly expressed on the surface of a variety of immune cells. Due to the additional carbohydrate chains, the molecular weight of CD37 fluctuates between approximately 40 and 64 kDa, CD37 is predominantly expressed in B-cells and is critical in regulating B-cell survival, shaping immune responses and immune evasion. CD37 plays a role in immune regulation and tumor suppression, and its immunotherapeutic targeting may provide additional benefit to patients with refractory B-cell malignancies. (PMID: 33333768) In acute myeloid leukemia (AML), CD37 has also been shown to be expressed on AML cells but not on normal hematopoietic stem cells (HSCs), which provides new therapeutic directions for the development of antibody-drug couplings (ADCs) targeting CD37.		
References	1. Bobrowicz M, et al. Int J Mol Sci. 2020;21(24):9531. 2. Caulier B, et al. Cell Rep Med. 2024;5(6):101572. 3. van Deventer S, et al. Biophys J. 2024;123(15):2301-2311. 4. Xu-Monette ZY, et al. Blood. 2016;128(26):3083-3100.		
Synonyms	CD37, Leukocyte antigen CD37		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Mouse CD37 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) conditions and stained using Coomassie blue.