

For Research Use Only

Recombinant Mouse Galanin protein (rFc Tag)



Catalog Number: Eg2758


Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Mouse Galanin protein Ala20-Ser124 (Accession# P47212) with a rabbit IgG Fc tag at the C-terminus. GeneID: 14419 Accession: P47212 Predicted Molecular Mass: 37.5 kDa SDS-PAGE: 36-42 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Galanin is a neuroendocrine peptide that is widely expressed in the central and peripheral nervous systems and also the gastrointestinal tract, pancreas, adrenal gland and urogenital tract. The encoded protein is a precursor that is proteolytically processed to generate two mature peptides: galanin and galanin message-associated peptide (GMAP). Galanin has diverse physiological functions including nociception, feeding and energy homeostasis, osmotic regulation, and water balance. GMAP has been demonstrated to possess antifungal activity and is hypothesized to be part of the innate immune system		
References	1. W E Schmidt, et al. (1991). Proc Natl Acad Sci U S A. 88(24):11435-9. 2. Sipin Zhu, et al. (2022). Front Endocrinol (Lausanne). 13:1019943. 3. J Šipková, et al. (2017). Physiol Res. 66(5):729-740.		
Synonyms	Gal, Galanin, Galanin message-associated peptide, Galanin peptides, Galn		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Mouse Galanin was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.