

For Research Use Only


Recombinant Mouse Surfactant protein D protein (rFc Tag)



Catalog Number: Eg2924

Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	<p>Purity: >90 %, SDS-PAGE</p> <p>Endotoxin Level: <0.1 EU/μg protein, LAL method</p> <p>Source: HEK293-derived Mouse Surfactant protein D protein Ala20-Phe374 (Accession# P50404) with a rabbit IgG Fc tag at the C-terminus.</p> <p>GeneID: 20390</p> <p>Accession: P50404</p> <p>Predicted Molecular Mass: 61.7 kDa</p> <p>SDS-PAGE: 67-75 kDa, reducing (R) conditions</p> <p>Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.</p>		
Biological Activity	Not tested		
Storage and Shipping	<p>Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.</p> <ul style="list-style-type: none">• Until expiry date, -20°C to -80°C as lyophilized proteins.• 3 months, -20°C to -80°C under sterile conditions after reconstitution. <p>Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.</p>		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	<p>Surfactant Protein D (SP-D) is a pattern recognition molecule that belongs to the collectin family expressed in multiple human organ systems, including alveolar type II cells and Clara cells. The surfactant protein D-encoding gene (SFTPD) is located at the genomic position 10q22.2-23.1. SP-D plays an important role in the innate immune system by binding bacteria, viruses, fungi, and parasites for clearance via opsonization in phagocytes, as well as aiding in the removal of allergens.</p>		
References	<ol style="list-style-type: none">1. Sorensen GL. (2018) Front Med. 5:18.2. Sorensen GL, et al. (2014) J Pediatr. 165(4):683-689.3. Wright JR. (2005) Nat Rev Immunol. 5(1):58-68.		
Synonyms	Lung surfactant protein D, PSP-D, Pulmonary surfactant-associated protein D, Sftp4, Sftpd		
For technical support and original validation data for this product please contact		This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.	
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)		E: proteintech@ptglab.com W: ptglab.com	

Selected Validation Data



Purity of Recombinant Mouse Surfactant protein D was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.