

For Research Use Only

Recombinant Human NBL1 protein (rFc Tag)(HPLC verified)



Catalog Number: Eg2940

Basic Information

Species:
Human

Purity:
>90 %, SDS-PAGE
 >90 %, SEC-HPLC

Tag:
rFc Tag

Technical Specifications

Purity:
>90 %, SDS-PAGE
 >90 %, SEC-HPLC

Endotoxin Level:
<0.1 EU/μg protein, LAL method

Source:
HEK293-derived Human NBL1 protein Ala17-Asp181 (Accession# P41271-1) with a rabbit IgG Fc tag at the C-terminus.

GeneID:
4681

Accession:
P41271-1

Predicted Molecular Mass:
43.7 kDa

SDS-PAGE:
45-55 kDa, reducing (R) conditions

Formulation:
Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

NBL1 (neuroblastoma suppression of tumorigenicity 1), also named as DAN, is the founding member of the DAN subfamily. It was originally identified as a putative tumor suppressor gene in a transformed fibroblast rat model, whereby it was significantly down-regulated in v-src-transformed rat fibroblasts. The MW of this protein is 19 kDa, and this antibody specially recognises the 19 kDa protein.

References

1. Tsai CH. et al. (2023) Biomolecules. 13(11):1570.
2. Ladenstein R. et al. (2018) Lancet Oncol. 19(12):1617-1629.
3. Lu L. et al. (2022) Oncogenesis. 11(1):21.
4. Du J. et al. (2023) Thorac Cancer. 14(30):3032-3041.
5. Wang A. et al. (2022) Parasitol Res. 121(5):1369-1378.

Synonyms


D1S1733E, DAN, DAND1, NB, Neuroblastoma suppressor of tumorigenicity 1

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Human NBL1 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Human NBL1 was greater than 90% as determined by SEC-HPLC.