

For Research Use Only

Recombinant Mouse IL-28B/ IFNL3 protein (rFc Tag)

Catalog Number: Eg3146




Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/µg protein, LAL method Source: HEK293-derived Human IL-28B protein Asp20-Val193 (Accession# Q8CGK6) with a rabbit IgG Fc tag at the C-terminus. GeneID: 338374 Accession: Q8CGK6 Predicted Molecular Mass: 45.6 kDa SDS-PAGE: 46-55 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	IL-28B, also named as IFNL3, is a cytokine distantly related to type I interferons and the IL-10 family. IL-28B, interleukin 28A (IL-28A), and interleukin 29 (IL-29) are three closely related cytokine genes that form a cytokine gene cluster on a chromosomal region mapped to 19q13. Expression of the cytokines encoded by the three genes can be induced by viral infection. All three cytokines have been shown to interact with a heterodimeric class II cytokine receptor that consists of interleukin 10 receptor, beta (IL-10RB) and interleukin 28 receptor, alpha (IL-28RA). IL-28B was identified as a key regulator of B- and T-cell vaccine responses against influenza.		
References	1. Sheppard P. et al. (2003). Nat Immunol. 4: 63-68. 2. Katrin Witte. et al. (2010). Cytokine Growth Factor Rev. 21(4):237-251. 3. Christabel Kelly. et al. (2011). Gut. 60(9):1284-1293. 4. Adrian Egli. et al. (2014). PLoS Pathog. 10(12):e1004556.		
Synonyms	Ifnl3, IFN-lambda-3, IL28, IL28b, IL-28B		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Mouse IL-28B was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) conditions and stained using Coomassie blue.