

For Research Use Only

Recombinant Mouse LBP protein (rFc Tag)



Catalog Number: Eg3163


Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Mouse LBP protein Gly25-Val481 (Accession# Q61805) with a rabbit IgG Fc tag at the C-terminus. GeneID: 16803 Accession: Q61805 Predicted Molecular Mass: 76.7 kDa SDS-PAGE: 80-90 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	LBP (Lipopolysaccharide-binding protein) is a 65 kDa type I acute response protein found in high concentrations in the blood and known to be of hepatic origin with antioxidant properties. It acts as an antioxidant to control lipid homeostasis and protects against oxidative stress by coupling with redox signaling and lipid metabolism. (PMID: 38615060) LBP inhibits monoamine biosynthesis by acting as an endogenous inhibitor of dopamine-β-hydroxylase (DBH) and aromatic-L-amino acid-decarboxylase (DDC), and targeting LBP may have therapeutic potential for some patients with depression.		
References	1. Zhang Q, et al. Nat Commun. 2024; 3:15(1):3213. 2. Fang M, et al. Immunity. 2023;56(3):620-634.e11. 3. Milbank E, et al. Pharmacol Res. 2023;187:106562.		
Synonyms	Lipopolysaccharide-binding protein		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Mouse LBP was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.