

For Research Use Only

Recombinant Rabbit CD8B protein (His Tag)



Catalog Number: Eg3262


Basic Information	Species: Rabbit	Purity: >90 %, SDS-PAGE	Tag: His Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Rabbit CD8B protein Ser19-Cys164 (Accession# G1SX39) with a His tag at the C-terminus. GeneID: 100355767 Accession: G1SX39 Predicted Molecular Mass: 17.2 kDa SDS-PAGE: 25-32 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	CD8 is a transmembrane glycoprotein predominantly expressed on cytotoxic T cells and can also be found on natural killer cells, cortical thymocytes, and dendritic cells. CD8b (T-cell surface glycoprotein CD8 beta chain) is an integral membrane glycoprotein that forms disulfide-linked heterodimers with CD8a (CD8 alpha chain). The CD8 alpha/beta heterodimer is the predominant CD8 complex expressed on the cell surface .		
References	1. Hennecke S. et al. (1993). J Biol Chem. 268(35):26607-12. 2. Xie Q, et al. (2021). Acta Pharm Sin B. 11(6):1365-1378. 3. Wheeler CJ. et al. (1992). Nature. 357(6375):247-9. 4. Terry LA. et al. (1990). Tissue Antigens. 35(2):82-91.		
Synonyms	CD8 subunit beta		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Rabbit CD8B was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.