

For Research Use Only

Recombinant Rat MMP-8 protein (mFc Tag)



Catalog Number: Eg3277

Basic Information

Species:
Rat

Purity:
>90 %, SDS-PAGE

Tag:
mFc Tag

Technical Specifications

Purity:

>90 %, SDS-PAGE

Endotoxin Level:

<0.1 EU/μg protein, LAL method

Source:

HEK293-derived Rat MMP-8 protein Leu21-Pro466 (Accession# O88766) with a mouse IgG Fc tag at the C-terminus.

GeneID:

63849

Accession:

O88766

Predicted Molecular Mass:

77.5 kDa

SDS-PAGE:

70-90 kDa, reducing (R) conditions

Formulation:

Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

Matrix metalloproteinase-8 (MMP-8) is a zinc-dependent endopeptidase that plays a crucial role in various biological processes, including tissue repair, inflammation, and neurodegenerative diseases. In rats, MMP-8 has been extensively studied, and several key findings have emerged. Recent evidence suggests that reactive astrocytes play an important role in neuroinflammation and neurodegenerative diseases. A study investigated whether an MMP-8 inhibitor (M8I) could control neuroinflammation in lipoteichoic acid (LTA)-stimulated rat primary astrocytes. The results showed that M8I significantly inhibited LTA-induced expression of inflammatory molecules such as iNOS, TNF-α, IL-1β, IL-6, and TLR-2. Additionally, M8I increased the expression of phase II antioxidant enzymes such as hemeoxygenase-1, NQO1, catalase, and MnSOD by modulating the Nrf2/ARE signaling pathway. These findings suggest the therapeutic potential of an MMP-8 inhibitor in neuroinflammatory disorders associated with astrocyte reactivity.

References

1. Yu J, et al. (2019) Arch Oral Biol. 97:238-244.
2. Tsubota M, et al. (2002) J Dent Res. 81(10):673-8.
3. Lee EJ, et al. (2018) J Neuroinflammation. 15(1):326.
4. Cederqvist K, et al. (2006) Pediatr Res. 60(4):395-400.

Synonyms


MMP8, EC:3.4.24.34, Matrix metalloproteinase-8, Neutrophil collagenase

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Rat MMP-8 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.