

For Research Use Only


Recombinant Rat CXCL7 protein (mFc Tag)



Catalog Number: Eg3404

Basic Information	Species: Rat	Purity: >90 %, SDS-PAGE	Tag: mFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Rat CXCL7 protein Ile46-Ile107 (Accession# A6KKE5) with a mouse IgG Fc tag at the N-terminus. GeneID: 246358 Accession: A6KKE5 Predicted Molecular Mass: 34.5 kDa SDS-PAGE: 34-37 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	CXCL7 is a platelet-derived growth factor that belongs to the CXC chemokine family, functioning as a potent chemo-attractant and activator of neutrophils through binding to its receptor CXCR2. CXCL7 has been demonstrated to participate in a variety of cellular processes, such as DNA synthesis, glycolysis, mitosis, intracellular cAMP accumulation, prostaglandin E2 secretion, as well as the synthesis of hyaluronic acid and plasminogen activator. In cancer cells, CXCL7 exerts promotion of cancer development in a variety of different cancer types, including renal cell carcinoma, lung cancer, and breast cancer.		
References	1. Yu M, et al. (2010) J Oncol. 2010:939407. 2. Guo Q, et al. (2017) Oncol Rep. 37(2):1114-1122. 3. Kinouchi T, et al. (2017) Cancer Sci. 108(12):2495-2502. 4. Wang YH, et al. (2021) Cell Death Dis. 12(12):1090.		
Synonyms	Chemokine (C-X-C motif) ligand 7, Ppbb		
For technical support and original validation data for this product please contact			
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)		E: proteintech@ptglab.com W: ptglab.com	
This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.			

Selected Validation Data



Purity of Recombinant Rat CXCL7 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.