

For Research Use Only

Recombinant Mouse CD314/NKG2D protein (rFc Tag)



Catalog Number: Eg3829

Basic Information

Species:
Mouse

Purity:
>90 %, SDS-PAGE

Tag:
rFc Tag

Technical Specifications

Purity:

>90 %, SDS-PAGE

Endotoxin Level:

<0.1 EU/μg protein, LAL method

Source:

HEK293-derived Mouse CD314 protein Phe90-Val232 (Accession# O54709-1) with a rabbit IgG Fc tag at the N-terminus.

GeneID:

27007

Accession:

O54709-1

Predicted Molecular Mass:

43.6 kDa

SDS-PAGE:

45-55 kDa, reducing (R) conditions

Formulation:

Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

Natural killer group 2D (NKG2D) (CD314) is an activating receptor on NK cells, CD8+ T cells, NKT cells and γδ+ T cells. Interleukin (IL)-15-induced up-regulation of NKG2D on intraepithelial cytotoxic T cells (IELs) and NKG2D ligands on enterocytes is a hallmark in CD and is responsible for epithelial cell destruction and subsequent villous atrophy.

References

1. Guy-Grand D, et al. Eur J Immunol. 1996;26:2248-2256.
2. Abadie V, et al. Semin Immunopathol. 2012;34:551-566.

Synonyms

CD314, CD314/ NKG2D, Klrk1, Killer cell lectin-like receptor subfamily K member 1, NK cell receptor D


For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Mouse CD314 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.