

For Research Use Only

Recombinant Human SAA2 protein (rFc Tag)



Catalog Number: Eg3837


Basic Information	Species: Human	Purity: >85 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >85 %, SDS-PAGE Endotoxin Level: <0.1 EU/µg protein, LAL method Source: HEK293-derived Human SAA2 protein Arg19-Tyr122 (Accession# P0DJ19-1) with a rabbit IgG Fc tag at the N-terminus. GeneID: 6289 Accession: P0DJ19-1 Predicted Molecular Mass: 38.8 kDa SDS-PAGE: 37-43 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Human SAA proteins are a group of apolipoproteins mainly found in the high density lipoprotein (HDL) portion of plasma. The serum amyloid A (SAA) protein is an acute phase apolipoprotein reactant produced mainly by hepatocytes and under the regulation of inflammatory cytokines. Reactive, secondary amyloidosis is characterized by the extracellular accumulation in various tissues of the SAA2 protein.		
References	1. Lundmark K, et al. (2002) Proc Natl Acad Sci U S A. 14;99(10):6979-84. 2. Baba S, et al. (1992) Biochim Biophys Acta. 10;1180(2):195-200.		
Synonyms	Amyloid A2 protein, serum amyloid A2, Serum amyloid A-2 protein		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Human SAA2 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.