

For Research Use Only


Recombinant Human Entactin protein (rFc Tag)



Catalog Number: Eg3912

Basic Information	Species: Human	Purity: >90 %, SDS-PAGE	Tag: N-rFc
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Human Entactin protein Asp385-Asn670 (Accession# P14543-1) with a rabbit IgG Fc tag at the N-terminus GeneID: 4811 Accession: P14543-1 Predicted Molecular Mass: 58.8 kDa SDS-PAGE: 55-65 kDa, reducing (R) condition Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Nidogen 1 (Entactin), first identified from a basement membrane-secreting cell line and the murine EHS tumor, is together with perlecan, laminin, and collagen IV, a ubiquitous component of basement membranes. Nidogen-1 is made by mesenchymal cells and is incorporated into the epithelial BM. It is a sulfated glycoprotein of 150 kDa with three globular domains, G1, G2 and G3, and contains binding sites for other ECM molecules.		
References	1. Murshed M, et al. (2000). Mol Cell Biol. Sep;20(18):7007-7012. 2. Pujuguet P, et al. (2000). J Cell Sci. Mar;113 (Pt 5)(Pt 5):849-858. 3. Vasudevan A, et al. (2010). Hippocampus. May;20(5):608-20.		
Synonyms			
For technical support and original validation data for this product please contact			
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)		E: proteintech@ptglab.com W: ptglab.com	
This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.			

Selected Validation Data



Purity of Recombinant Human Entactin was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.