

For Research Use Only

Recombinant Mouse VISTA protein (mFc Tag)(HPLC verified)



Catalog Number: Eg4197

Basic Information

Species:
Mouse

Purity:
>90 %, SDS-PAGE
 >90%, SEC-HPLC

Tag:
mFc Tag

Technical Specifications

Purity:
>90 %, SDS-PAGE
 >90%, SEC-HPLC

Endotoxin Level:
<0.1 EU/μg protein, LAL method

Source:
HEK293-derived Human VISTA protein Phe33-Ala191 (Accession# Q9D659) with a mouse IgG Fc tag at the C-terminus.

GeneID:
74048

Accession:
Q9D659

Predicted Molecular Mass:
44.4 kDa

SDS-PAGE:
50-65 kDa, reducing (R) conditions

Formulation:
Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

VISTA, also known as PD-1H, Dies1, and Gi24, is encoded by the Vsir gene in mice. VISTA is a type I transmembrane protein that acts as an immune-regulatory checkpoint on hematopoietic cells and shares homology with both the B7 family ligand PD-L1 and the CD28 family receptor PD-1. In mice, VISTA is expressed predominantly in hematopoietic tissues, such as the spleen, thymus, lymph node, and bone marrow (BM). VISTA plays a key role in regulating immune system responses and maintaining homeostasis.

References

- 1.Yuan L, et al. (2021) Trends Immunol. 42(3):209-227.
- 2.Tagliamento M, et al. (2021) Immunotargets Ther. 10:185-200.
- 3.Martin AS, et al. (2023) Front Immunol. 14:1086102.
- 4.Shekari N, et al. (2023) Cancer Cell Int. 23(1):265.
- 5.Zhang RJ, et al. (2024) Exp Mol Med. 56(11):2348-2356.

Synonyms


Dies1, PD-1H, Platelet receptor Gi24, V-set domain-containing immunoregulatory receptor, V-set immunoregulatory receptor

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Mouse VISTA was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Mouse VISTA was greater than 90% as determined by SEC-HPLC.