

For Research Use Only

Recombinant Canine SCF protein (rFc Tag)



Catalog Number: Eg4548


Basic Information	Species: Canine	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/µg protein, LAL method Source: HEK293-derived Canine SCF protein Lys26-Gln215 (Accession# Q06220) with a rabbit IgG Fc tag at the C-terminus. GeneID: 403507 Accession: Q06220 Predicted Molecular Mass: 47.3 kDa SDS-PAGE: 50-65 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Stem cell factor (SCF), also known as mast cell growth factor, steel factor and Kit ligand, is a haematopoietic growth factor. SCF is a critical protein with key roles in the cell such as haematopoiesis, gametogenesis and melanogenesis. SCF activates multiple signalling pathways by binding to the c-kit receptor. SCF and its receptor, c-kit ligand (KL), are particularly upregulated in human malignancies, including gastrointestinal stromal tumours (GISTs), breast cancer, haematopoietic malignancies, myeloid leukaemia and gliomas. The SCF/c-KIT system also plays an important role in cell fate decision and its activity is known to control cell proliferation, differentiation, survival and apoptosis.		
References	1.Vazquez-Mellado MJ. et al. (2017) Vitam Horm. 105:273-296. 2.Figueira MI. et al. (2017) Cancer Lett. 405:10-21. 3.Yasuda A. et al. (2006) Mol Cancer. 5:46.		
Synonyms	KITLG, c-Kit ligand, Kit ligand, Mast cell growth factor, MGF		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Canine SCF was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.