

For Research Use Only

Recombinant Human HSP27 protein (rFc Tag)



Catalog Number: Eg4620


Basic Information	Species: Human	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Human HSP27 protein Ala70-Thr184 (Accession# P04792) with a rabbit IgG Fc tag at the C-terminus. GeneID: 3315 Accession: P04792 Predicted Molecular Mass: 39.0 kDa SDS-PAGE: 40-55 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Heat shock protein 27 (HSP27), also known as HSPB1, belongs to the small heat shock protein family which is induced in response to environmental challenges or/and developmental transitions. It is also an anti-apoptotic protein that plays crucial roles in tumorigenesis and cell survival and is reported to be an independent prognosis marker for cancer. Recently HSP27 has been found to be a valuable marker for melanoma.		
References	1. Kijima K. et al. (2005). J Hum Genet. 50(9):473-476. 2. Foster CS. et al. (2009). Br J Cancer. 101(7):1137-44. 3. Okuno M. et al. (2016). Int J Mol Sci. 17(1):137. 4. Seo MK. et al. (2020). Pharmacoecon Open. 4(3):529-539.		
Synonyms	HSPB1, CMT2F, Estrogen-regulated 24 kDa protein, HMN2B, Hsp25		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Human HSP27 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.