

For Research Use Only

Recombinant Mouse Ferritin heavy chain/FTH1 protein (rFc Tag)



Catalog Number: Eg4765


Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/μg protein, LAL method Source: HEK293-derived Mouse Ferritin heavy chain protein Met1-Ser182 (Accession# P09528) with a rabbit IgG Fc tag at the C-terminus. GeneID: 14319 Accession: P09528 Predicted Molecular Mass: 47.3 kDa SDS-PAGE: 50-60 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Ferritin heavy chain (FTH1, ferritin, heavy polypeptide 1), as a 21 kDa subunit of the ferritin complex, exhibits ferroxidase activity, which serves an essential role in catalyzing the conversion of the ferrous ions (Fe²⁺) to the ferric form (Fe³⁺). FTH1 regulates angiogenesis during inflammation and malignancy by binding to high molecular weight kininogen (HK) and blocking its anti-angiogenic activity on endothelial cells in vitro and in vivo. FTH1 has been shown to protect proximal tube epithelial cells and kidneys against the activity of free iron in reactive oxygen species generation. FTH1 plays an important role in the maintenance of the cellular iron balance in ferroptosis.		
References	1. Mimi C Sammarco, et al. (2008). J Biol Chem. 283(8):4578-87. 2. Maddalena Di Sanzo, et al. (2020). Cells. 9(12):2554.		
Synonyms	EC:1.16.3.1, Ferritin H subunit, Ferritin heavy chain, Ferritin heavy chain, N-terminally processed, Fth		

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Purity of Recombinant Mouse Ferritin heavy chain was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) conditions and stained using Coomassie blue.