

For Research Use Only

FITC Plus-conjugated CD31 Polyclonal antibody

Catalog Number: FITC-28083



Basic Information

Catalog Number: FITC-28083	GenBank Accession Number: NM_001032378	Purification Method: Antigen affinity purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 18613	Recommended Dilutions: IF-P: 1:50-1:500
Source: Rabbit	UNIPROT ID: Q08481	Excitation/Emission maxima wavelengths: 495 nm / 524 nm
Isotype: IgG	Full Name: platelet/endothelial cell adhesion molecule 1	
Immunogen Catalog Number: AG27828	Calculated MW: 81 kDa	
	Observed MW: 110-130 kDa	

Applications

Tested Applications: IF-P	Positive Controls: IF-P: mouse brain tissue,
Species Specificity: mouse	

Background Information

Platelet endothelial cell adhesion molecule-1 (PECAM-1, CD31) is a member of the immunoglobulin gene superfamily of cell adhesion molecules. CD31 is a transmembrane glycoprotein that is highly expressed on the surface of the endothelium, making up a large portion of its intracellular junctions. PECAM-1 is also present on the surface of hematopoietic cells and immune cells including platelets, monocytes, neutrophils, natural killer cells, megakaryocytes and some types of T-cell (PMID: 9011572). As well as its role in cell-cell adhesion, PECAM-1 functions as a signaling receptor, and is involved in important physiological events such as nitric oxide production, regulation of T-cell immunity and tolerance, leukocyte transendothelial migration and inflammation and angiogenesis (PMID: 21183735; 20978210; 17872453; 20634489).


Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse brain tissue using FITC Plus CD31 antibody (FITC-28083) at dilution of 1:200. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).