

For Research Use Only

FITC Anti-Mouse CD115 (c-fms) (AFS98)



Catalog Number: FITC-65095

Basic Information

Catalog Number:

FITC-65095

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a

GenBank Accession Number:

BC036343

GeneID (NCBI):

12978

Full Name:

colony stimulating factor 1 receptor

Purification Method:

Affinity purification

CloneNo.:

AFS98

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ Balb/c mouse peritoneal macrophages were surface stained with 0.5 ug FITC Anti-Mouse CD115 (c-fms) (FITC-65095, Clone: AFS98) (green) or 0.5 ug FITC-Rat IgG2a Isotype Control (black). Samples were not fixed.