

For Research Use Only

FITC Plus Anti-Mouse CD41 (MWReg30)



Catalog Number: FITC-65189

Basic Information

Catalog Number: FITC-65189	GenBank Accession Number: BC120493	Purification Method: Affinity purification
Size: 100ug, 0.5 mg/ml	GeneID (NCBI): 16399	CloneNo.: MWReg30
Source: Rat	Full Name: integrin alpha 2b	Excitation/Emission maxima wavelengths: 495 nm / 524 nm
Isotype: IgG1, kappa		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse peripheral blood platelets were surface stained with 0.5 ug FITC Plus Anti-Mouse CD41 (FITC-65189, Clone: MWReg30) (red), or 0.5 ug Isotype Control. Cells were not fixed.



1X10⁶ mouse peripheral blood platelets were surface stained with 0.5 ug FITC Plus Anti-Mouse CD41 (FITC-65189, Clone: MWReg30) (red), or 0.5 ug Isotype Control. Cells were not fixed.